# BioRoot™ RCS

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### Josette Camilleri

B.Ch.D., M.Phil., Ph.D., FICD, FADM, FIMMM, FHEA (UK)  
School of Dentistry,  
Institute of Clinical Sciences  
College of Medical and Dental Sciences  
The University of Birmingham  
Birmingham  
U.K.

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## Literatur

1. Schilder H. Filling root canals in three dimensions. *Dent Clin North Am.* 1967.
2. Grossman LI. *Endodontic Practice.* Philadelphia: Lea & Febiger.1978
3. Sundqvist G, Figdor D. Endodontic treatment of apical periodontitis. In: Ørstavik D, Pitt Ford TR, eds. *Essential Endodontology. Prevention and Treatment of Apical Periodontitis.* Oxford: Blackwell, 1998.
4. Schembri M, Peplow G, Camilleri J. Analyses of heavy metals in mineral trioxide aggregate and Portland cement. *J Endod.* 2010;36(7):1210-5.
5. Chang SW, Shon WJ, Lee W, Kum KY, Baek SH, Bae KS. Analysis of heavy metal contents in gray and white MTA and 2 kinds of Portland cement: a preliminary study. *Oral Surg Oral Med Oral Pathol Oral Radiol Endod.* 2010;109(4):642-6.
6. Monteiro Bramante C, Demarchi AC, de Moraes IG, Bernadineli N, Garcia RB, Spångberg LS, Duarte MA. Presence of arsenic in different types of MTA and white and gray Portland cement. *Oral Surg Oral Med Oral Pathol Oral Radiol Endod.* 2008;106(6):909-13.
7. Camilleri J. Characterization and hydration kinetics of tricalcium silicate cement for use as a dental biomaterial. *Dent Mater.* 2011;27(8):836-44.
8. Xuereb M, Vella P, Damidot D, Sammut CV, Camilleri J. In situ assessment of the setting of tricalcium silicate-based sealers using a dentin pressure model. *J Endod.* 2015;41(1):111-24.
9. Demirkaya K, Can Demirdöğen B, Öncel Torun Z, Erdem O, Çetinkaya S, Akay C (2016). In vivo evaluation of the effects of hydraulic calcium silicate dental cements on plasma and liver aluminium levels in rats. *Eur J Oral Sci.* 124(1):75-81.
10. Demirkaya K, Demirdöğen BC, Torun ZÖ, Erdem O, Çırak E, Tunca YM (2016). Brain aluminium accumulation and oxidative stress in the presence of calcium silicate dental cements. *Hum Exp Toxicol.* pii: 0960327116679713.
11. Forbes WF, Gentleman JF. Risk factors, causality, and policy initiatives: the case of aluminum and mental impairment. *Exp Gerontol* 1998;33:141–54.
12. Arias-Moliz MT, Camilleri J. The effect of the final irrigant on the antimicrobial activity of root canal sealers. *J Dent.* 2016;52:30-6.
13. Viapiana R, Moizadeh AT, Camilleri L, Wesselink PR, Tanomaru Filho M, Camilleri J. Porosity and sealing ability of root fillings with gutta-percha and BioRoot™ RCS or AH Plus sealers. Evaluation by three ex vivo methods. *Int Endod J.* 2016;49(8):774-82.
14. Collado-González M, García-Bernal D, Oñate-Sánchez RE, Ortolani-Seltenerich PS, Lozano A, Forner L, Elena C, Rodríguez-Lozano FJ. Biocompatibility of three new calcium silicate-based endodontic sealers on human periodontal ligament stem cells. *Int Endod J.* 2016 Sep 26. doi: 10.1111/iej.12703. [Epub ahead of print].
15. Poggio C, Riva P, Chiesa M, Colombo M, Pietrocola G. Comparative cytotoxicity evaluation of eight root canal sealers. *J Clin Exp Dent.* 2017;1;9(4):e574-e578.
16. Camps J, Jeanneau C, El Ayachi I, Laurent P, About I. Bioactivity of a Calcium Silicate- based Endodontic Cement (BioRoot™ RCS): Interactions with Human Periodontal Ligament Cells In Vitro. *J Endod.* 2015;41(9):1469-73.
17. Camilleri J. Sealers and warm gutta-percha obturation techniques. *J Endod.* 2015;41(1):72-8.
18. Jeong JW, DeGraft-Johnson A, Dorn SO, Di Fiore PM. Dentinal Tubule Penetration of a Calcium Silicate-based Root Canal Sealer with Different Obturation Methods. *J Endod.* 2017;43(4):633-637.
19. McMichael GE, Primus CM, Opperman LA. Dentinal Tubule Penetration of Tricalcium Silicate Sealers. *J Endod.* 2016;42(4):632-6.
20. Alshehri M, Alamri HM, Alshwaimi E, Kujan O. Micro-computed tomographic assessment of quality of obturation in the apical third with continuous wave vertical compaction and single match taper sized cone obturation techniques. *Scanning.* 2016;38(4):352-6.
21. Donnermeyer D, Bunne C, Schäfer E, Dammaschke T. Retreatability of three calcium silicate-containing sealers and one epoxy resin-based root canal sealer with four different root canal instruments. *Clin Oral Investig.* 2017 Jun 22. doi: 10.1007/ s00784-017-2156-5. [Epub ahead of print]

**Septodont GmbH**  
Felix-Wankel-Str. 9  
53859 Niederkassel, Deutschland  
Tel.: +49 (0) 228 971 26-0  
Fax: +49 (0) 228 971 26-66  
E-Mail: [info@septodont.de](mailto:info@septodont.de)  
[www.septodont.com](http://www.septodont.com)